

Max3™ Narrow Breast Board

The versatile positioning device for upper body treatments.

BACK

CT friendly and IMRT ready, the NEW Carbon Fiber Max3 Narrow Breast Board is 16" at the widest point. Minimal loose parts and ergonomic design help to maximize function, ease of use and patient comfort.

Carbon Fiber wedge angles are now connected and stored in the breast board base. These wedge angles are designed to distribute weight evenly across the whole board, minimizing any bowing or flexing while adding to the strength of the board. The Max3 is surprisingly lightweight and easy to transport.

Take a look at the new features the Max3 has to offer:

The Breast Board offers many options and features.


- Accommodates Silverman Head supports.
- Tapered design at the lower back for comfortable transition.
- A built in adjustable locking Backstop cushion to promote reproducibility and comfort.
- Secure Velcro closure locks board to base for easy transport and storage.
- Accommodates Bionix reusable and disposable head only thermoplastic masks.
- Accommodates the use of Butterfly Boards.
- Elevation wedges (7°, 12°, 17°, 23°). Custom wedges can be added.
- Locks down to most couches.

Max 3 includes head cup, hand post, thermoplastic handles, angles and Backstop.

Service

All Bionix Radiation Therapy products are fully guaranteed and warranted to be free from defects in workmanship for a period of three years from date of purchase.

*To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.*



Optional calibrated arm and wrist cups securely lock in place.



BACK

United States Patent Number 406,345

