

SuProne Plus™ System

The durable positioning device for head and neck treatments.

The NEW SuProne Plus Head & Neck System allows secure, accurate and versatile positioning. This new system will allow you to position prone, supine, lateral and flat with maximum patient comfort. This cost efficient solution to all your head and neck treatments.

SuProne Plus Head & Neck

- Fits both disposable or reusable frames
- Timo/Silverman or MoldCare® attachment
- Locks down to Exact™ Couch
- Prone or supine treatment
- Lateral attachment
- 45° attachment
- Prone attachment
- May be used for flat positioning
- Compatible with both 9" x 10" and 10" x 10" reusable frames
- Includes a comfort cushion for patient in prone position

Service

All Bionix Radiation Therapy products are fully guaranteed and warranted to be free from defects in workmanship for a period of three years from date of purchase.

To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.

****See following page for attachments to maximize your treatment potentials***

RT-5039



SuProne Plus Head & Neck

RT-7039



SuProne Plus in Carbon Fiber



SuProne Plus™ System Attachments

Attachments for the durable positioning device.

The NEW SuProne Plus Head & Neck System allows secure, accurate and versatile positioning. This new system will allow you to position prone, supine, lateral and flat with maximum patient comfort. This cost efficient system allows one board with attachments to solve all your head and neck treatments.

NEW 45° Attachment

- Attaches to SuProne Plus Head & Neck Board (RT-5039 & RT-7039)
- May be positioned 5 degrees to 45 degrees
- Easy to locate and read degree location
- Stable base for maximum reproducibility
- Locks down to most couches
- Attaches to VersaBoard™

Prone Attachment (also available in acrylic - RT-5029)

- Attaches to SuProne Plus Head & Neck Board (RT-5039 & RT-7039)
- May be positioned -20 degrees or +20 degrees
- Includes comfort cushion for patient
- Easy to locate and read degree location
- Stable base for maximum reproducibility
- Locks down to most couches
- Attaches to VersaBoard™

Lateral Attachment

- Attaches to SuProne Plus™ Head and Neck Board (RT-5039 & RT-7039)
- Three positioning heights available (2cm, 5cm and 7cm)
- Comfortable lateral positioning
- May be used prone
- Locks down to most couches
- Attaches to VersaBoard™

Service

All Bionix Radiation Therapy products are fully guaranteed and warranted to be free from defects in workmanship for a period of three years from date of purchase.

To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.



SuProne Plus™ shown with 45° attachment



SuProne Plus™ shown with prone attachment



SuProne Plus™ shown with lateral attachment & SecureVac™

