

TorsoBoard SE

The cost effective positioning device for use in the simulation room.

The TorsoBoard SE is an acrylic version of our popular carbon fiber Max2™ TorsoBoard, but the TorsoBoard SE was designed to be used in the simulation room. The TorsoBoard SE has the same dimensions, positioning holes, and arm supports as the Max2™ so you can establish positioning protocol in the simulation room that can be used in the treatment room—without ever moving the device. And because it's made out of acrylic, you save money.

Versatile

The TorsoBoard SE is designed to be used for both lung and breast treatments and is compatible with low-melt thermoplastics used for immobilization.

Accuracy

Like all Bionix products, the TorsoBoard SE is manufactured using precision computer-aided machining to ensure you receive the highest quality product.

Other Features Include:

- SecureFit Wedges (5°, 10°, 15°, 20°, 25°)
- Accommodates Timo, Silverman, or Head Cup
- Skid Resistant Base
- Optional Multi-Purpose Arm Rest
- Optional Back Stop

Service

All Bionix Radiation Therapy products are fully guaranteed and warranted to be free from defects in workmanship for a period of three years from date of purchase.

*To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.*



RT-5045



The TorsoBoard SE is an acrylic duplicate of the Max2™



RT-5050

Optional Back Stop also available to help keep patient from slipping. Easily attaches to all Bionix Breast Boards with ease.

