

# Wedge Breast Board & Max Wedge Breast Board

Versatile positioning devices for breast and lung treatments.

## Bionix® Breast Boards

Breast boards are ideal for treatment of the breast or lungs. The use of low-melt thermoplastic material allows easy positioning and immobilization of the patient in either flat or inclined positions. Bionix used a wedge system for better stability and both boards include thermoplastic handles, set of 5 wedges, head cup, hand post and transition pad.

### RT-5000 Wedge Breast Board

- Made of Acrylic
- Uses Wedge System (5°, 10°, 15°, 20°, 25°)
- Includes hand post, transition pad, thermoplastic handles and head cup
- Attaches to most couches

### RT-5026 Max Wedge Breast Board

- Designed with carbon fiber panels to minimize attenuation
- Uses Wedge System 5-25
- Includes hand post, transition pad, thermoplastic handles and head cup
- Attaches to most couches

## Service

All Bionix Radiation Therapy products are fully guaranteed and warranted to be free from defects in workmanship for a period of three years from date of purchase.

To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at [www.BionixRT.com](http://www.BionixRT.com).



### RT-5050

Optional Back Stop also available to help keep patient from slipping. Easily attaches to all Bionix Breast Boards with ease.



The wedges distribute the weight evenly so that bowing does not occur. The board can be put in a flat position or inclined at five different degrees.

